

THE VALLEY HEALTH SYSTEM STUDY WITH THE TYRX™ ANTIBACTERIAL ENVELOPE*

High-Risk Patients Implanted with the TYRX Antibacterial Envelope are Significantly Less Likely to Develop Cardiac Implantable Electronic Device (CIED) Infection than Comparatively Low-Risk Cohorts¹

DESIGN

A retrospective, dual-cohort study was conducted in a large population of patients undergoing CIED procedures to determine the effect of the TYRX Antibacterial Envelope on CIED Infection rates, utilizing a novel scoring index to risk-stratify patients based on the specific combination of risk factors, rather than just the absolute number of risk factors.¹

6-MONTH CIED INFECTION RATE: STRATIFIED BY DEVICE AND PROCEDURE TYPE(%) ¹				
Device	Procedure Type			
	DeNovo	Generator	Upgrade	Other
CRT-D	4.0%	3.1%	3.9%	33.0%
ICD	1.5%	1.7%	18.2%	0.0%
PM†	0.8%	0.3%	0.0%	0.0%
CRT-P‡	0.0%	0.0%	0.0%	0.0%

Infection rates were higher in the Implantable Cardioverter Defibrillator (ICD) and Cardiac Resynchronization Therapy Defibrillator (CRT-D) groups, and in upgrade and early pocket re-exploration group.

† Pacemaker ‡ Cardiac Resynchronization Therapy Pacemaker

METHODS

Two cohorts of patients who underwent CIED procedures were identified: 1,651 patients before the introduction of the TYRX Envelope at the site (January 2007-October 2009) and 1,240 patients after the introduction of the TYRX Envelope (October 2009-September 2011), including 275 patients who received the TYRX Envelope. Using propensity-score matching, the 275 patients who received the TYRX Envelope were matched to 275 patients prior to the introduction of the TYRX Envelope.¹

RESULTS

Compared to patients who did not receive a TYRX Envelope, those who did receive the Envelope were more likely to have the following risk factors associated with an increased risk of an infection: early pocket re-exploration, male gender, diabetes, device upgrade, congestive heart failure (CHF), hypertension, GFR < 60ml/min. All these variables have been associated with CIED Infection.¹


The 6-month infection rate was significantly lower in patients who received a TYRX Envelope, compared to propensity score matched patients who did not (1.1% vs. 3.6%, p=0.048).¹

*Study performed utilizing the TYRX™ Non-Absorbable Antibacterial Envelope.

1. Mittal S et al. *Heart Rhythm*. 2014;11(4):595-601.

Medtronic TYRX, Inc. 1 Deer Park Drive, Suite G | Monmouth Junction, NJ | Medtronic CRHF/ Infection Control
TYRX Customer Service: 800.848.9300 Technical Support/Services: 800.505.4636 (Brady) and 800.723.4636 (Tachy) | www.TYRX.com

UC201503788bEN



A composite risk score was created by weight, adjusting for the seven risk factors: 3 groups emerged—low-risk (score 0-7: 1% infection), medium risk (score 8-14: 3.4% infection), and high risk (score ≥ 15: 11.1% infection).

CONCLUSIONS

The TYRX Envelope reduced infections by 79% and 100% in the medium and high-risk groups, respectively.¹

There were ~70% fewer infections in patients who received the TYRX Antibacterial Envelope, compared to those who did not, across all device types.¹

